

# Encapsulation Efficiency Analysis of RNA-LNPs Using RiboGreen Assay

2026

Protocol for fluorescence-based quantification of RNA in LNPs




# Introduction

The **Quant-iT™ RiboGreen™ RNA assay** (1) is a highly sensitive fluorescence-based method for RNA quantification. In the context of lipid nanoparticles (LNPs), it is commonly used to determine **encapsulation efficiency (EE%)** by distinguishing between **free (unencapsulated) RNA** and **total RNA** following nanoparticle disruption.

RiboGreen is a fluorescent dye that selectively binds to RNA, generating a signal proportional to RNA concentration. By measuring fluorescence before and after LNP lysis using a detergent (e.g., Triton X-100), it is possible to quantify:

- Free RNA (accessible before lysis)
- Total RNA (after LNP disruption)

Encapsulation efficiency is then calculated as the **proportion of RNA encapsulated within nanoparticles** (Fig. 1).



**Figure 1. Principle of encapsulation efficiency determination using the RiboGreen assay.** Before lysis, the RiboGreen dye selectively detects free (unencapsulated) RNA. Following disruption of LNPs with a detergent, the total RNA becomes accessible to the dye. Encapsulation efficiency is calculated from the difference between total and free RNA relative to the total RNA content. (2)

The RiboGreen assay is widely used in laboratories due to its **simplicity and compatibility with standard fluorescence microplate readers**, making it well suited for routine characterization of RNA-loaded nanoparticles. It offers **high sensitivity**, enabling detection of low RNA concentrations across a broad dynamic range. However, accurate quantification requires careful control of experimental conditions, as fluorescence signals may be influenced by the **sample environment, dilution, and measurement settings**.

## Key definitions & calculations

**Encapsulation efficiency (EE%)** reflects the proportion of RNA successfully encapsulated within nanoparticles and is a key indicator of formulation quality.

$$EE\% = \frac{[Total\ RNA - Free\ RNA]}{[Total\ RNA]} \times 100 \quad (Eq.1)$$

**Encapsulation yield (EY%)** represents the overall process efficiency. It may reflect potential losses occurring during formulation alone, or across both formulation and downstream processing steps, depending on the stage at which the measurement is performed.

$$EY\% = \frac{[Total\ RNA - Free\ RNA]}{[RNA_{input}]} \times 100 \quad (Eq.2)$$

In Eq.2, RNA input corresponds to the RNA concentration measured from the initial aqueous RNA solution used for formulation.

**RNA loading** describes the amount of RNA associated with the nanoparticle relative to lipid content and is relevant for dosing considerations. It can be expressed as weight percentage (wt%) or as µg RNA per mg of LNP.

$$RNA\ loading\ \% = \frac{Mass\ of\ RNA\ in\ LNPs}{Total\ mass\ of\ lipids\ in\ LNPs} \times 100 \quad (Eq.3)$$

## Materials

Table 1. List of equipment, reagents, and consumables required for EE% analysis using RiboGreen assay.

<b>Equipment</b>	<ul style="list-style-type: none"> <li>• Fluorescence microplate reader (e.g., TECAN)</li> <li>• Pipettes</li> </ul>
<b>Reagents</b>	<ul style="list-style-type: none"> <li>• Quant-iT™ RiboGreen™ RNA Reagent</li> <li>• TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 7.5)</li> <li>• Nuclease-free water</li> <li>• RNA standard</li> <li>• Detergent (e.g., Triton X-100, 10% v/v in nuclease-free water)</li> </ul>
<b>Consumables</b>	<ul style="list-style-type: none"> <li>• 96-well black microplate</li> <li>• Nuclease-free tubes and pipette tips</li> <li>• Plate cover (optional)</li> </ul>

# Protocol for EE% analysis of RNA-LNPs using RiboGreen

## 1. Preparation of reagents & solutions

### 1.1. Preparation of TE buffer

1. Prepare 1X TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 7.5) by diluting the provided 20X stock solution 20-fold with nuclease-free water.

**Note (!):** Ensure all solutions are **RNase-free** to prevent RNA degradation.

2. Use 1X TE buffer for:
  - Preparation of RiboGreen working solution
  - Dilution of RNA standards
  - Dilution of samples

### 1.2. Preparation of RiboGreen working solution (high-range assay)

3. Allow the RiboGreen reagent to equilibrate to room temperature before opening.
4. Prepare a 200-fold dilution in TE buffer (high-range assay, 20 ng/mL – 1 µg/mL RNA). Mix gently and protect the solution from light.
  - e.g., 50 µL RiboGreen + 9.95 mL TE buffer for ~100 samples

**Notes (!):**

- Prepare fresh and use within a few hours.
- Use **polypropylene tubes**; avoid glass containers (adsorption effects).
- Keep the reagent **protected from light** to prevent photodegradation.
- Based on prior internal studies, the **high-range RiboGreen assay** was selected for this protocol, as it provides optimal sensitivity and compatibility with RNA concentrations typically obtained in standard RNA-LNP formulations prepared using the TAMARA platform in our laboratory. (3)

### 1.3. Preparation of RNA stock solution

5. Prepare a 2 µg/mL RNA solution in TE buffer using RNase-free consumables.

**Notes (!):**

- We recommend preparing the standard curve using the **same RNA type and sequence as the one encapsulated in the LNPs**, rather than using the RNA

standard provided in the assay kit. This ensures more representative fluorescence behavior and improves the accuracy of quantification.

- The RNA standard solution should be prepared in the **same buffer conditions as the samples** to ensure accurate comparison.

## 1.4. Pre-dilution of RNA samples

6. Prior to plate setup, dilute RNA-LNP samples and RNA input solutions **1:20** in TE buffer.

### Notes (!):

- Following plate preparation (20  $\mu\text{L}$  pre-diluted sample + 80  $\mu\text{L}$  TE buffer + 100  $\mu\text{L}$  RiboGreen), the **final dilution factor reaches 1:200** relative to the initial sample concentration. This dilution must be taken into account during data analysis and when back-calculating RNA concentrations.
- The pre-dilution step is required to ensure that the final RNA concentrations fall within the **linear range of the high-range RiboGreen assay**. The pre-dilution factor may be **adjusted if necessary** depending on the RNA concentration of the formulation.

## 2. Plate setup & pipetting workflow

### 2.1. Deposition of RNA standards (high-range assay)

7. Dispense the **appropriate volumes** of RNA standard solution (2  $\mu\text{g}/\text{mL}$ ) into the designated wells according to **Table 2**. Prepare standards in **triplicate** across each row (e.g., A1–A3, B1–B3, etc.).

Table 2. Preparation of high-range RNA standard curve.

Well position	Volume of 2 $\mu\text{g}/\text{mL}$ RNA solution ( $\mu\text{L}$ )	Volume of TE buffer ( $\mu\text{L}$ )	Volume of RiboGreen ( $\mu\text{L}$ )	Final RNA concentration ( $\mu\text{g}/\text{mL}$ )
A1–A3	100	0	100	1
B1–B3	50	50	100	0.5
C1–C3	20	80	100	0.2
D1–D3	10	90	100	0.1
E1–E3	5	95	100	0.05
F1–F3	2	98	100	0.02
G1–G3	100	0	100	0

**Note (!):** The high-range assay is used in this protocol for standard curve preparation. For samples with lower RNA concentrations, the low-range assay may be used as an alternative. For low-range assay conditions, please refer to the Quant-iT™ RiboGreen™ RNA Reagent and Kit user guide (1).

## 2.2. Deposition of RNA-LNP samples

8. Add 20  $\mu\text{L}$  of each pre-diluted RNA-LNP sample (see Section 1.4) into the designated wells. Perform all measurements in **triplicate**.

### Notes (!):

- Samples are loaded such that a **final pre-dye volume of 100  $\mu\text{L}$  per well** is reached after TE addition (see Section 2.3). Under the conditions described in this protocol, samples are typically diluted to reach a **final RNA concentration of  $\sim 0.1 \mu\text{g}/\text{mL}$  in the assay well**, corresponding to the mid-range of the high-range standard curve for optimal accuracy. A typical approach for standard RNA-LNP formulations (e.g., SM-102, ALC-0315, LP-01) prepared using the Inside Therapeutics RNA-LNP Formulation Protocol with TAMARA (3) is to load **20  $\mu\text{L}$  of sample and add 80  $\mu\text{L}$  of TE buffer per well**.
- The **sample volume may be adjusted if necessary**, provided that the final pre-dye volume remains constant (100  $\mu\text{L}$ ). Higher dilution factors may help reduce the influence of interfering contaminants; however, very small sample volumes should be avoided due to reduced pipetting accuracy. (1)
- (Optional) For **encapsulation yield (EY%)** determination, include the **input RNA solution used during microfluidic mixing** in the plate. Load this solution into wells H1–H3 using the same conditions as samples (e.g., 20  $\mu\text{L}$  of pre-diluted input RNA solution).

## 2.3. Deposition of TE buffer

9. Add TE buffer to each well to reach a **final pre-dye volume of 100  $\mu\text{L}$** :
  - For **standards**: complete according to **Table 2**
  - For **samples**: add **80  $\mu\text{L}$  of TE buffer** to wells containing 20  $\mu\text{L}$  of sample

## 2.4. Deposition of RiboGreen reagent

10. Add 100  $\mu\text{L}$  of RiboGreen working solution to all wells (standards & samples). Ensure the plate is protected from light.

## 3. Fluorescence measurement – Free RNA

11. Place the microplate in the fluorescence plate reader. Set the instrument parameters as follows:

- **Temperature:** 25°C
  - **Shaking:** Brief mixing prior to reading (e.g., 30 seconds, linear shaking)
  - **Fluorescence intensity:**
    - **Excitation wavelength:** 485 nm
    - **Emission wavelength:** 530 nm
12. Measure fluorescence intensity for all wells.
- ❖ This first measurement corresponds to **free (unencapsulated) RNA** present outside the LNPs.

**Notes (!):**

- Set the instrument gain so that the wells with the highest RNA concentration produce fluorescence signals close to the maximum detection limit of the microplate reader.
- Ensure that all fluorescence values fall within the **linear range of the standard curve** for accurate quantification.
- Keep **measurement timing consistent** and **minimize exposure to light** to reduce photobleaching effects.

## 4. LNP disruption & Total RNA measurement

### 4.1. Addition of detergent (LNP lysis)

13. Add 10 µL of Triton X-100 (10% v/v) to each well (standards & samples).

14. Mix well to ensure homogeneous distribution of the detergent.

**Note (!):** Gentle pipette mixing during Triton X-100 addition can improve homogenization. However, avoid vigorous mixing to prevent the formation of bubbles, which may interfere with fluorescence measurements.

### 4.2. Incubation

15. Place the plate in the microplate reader and incubate for **25 minutes at 25°C with regular shaking (30 s shaking followed by 4 min 30 s rest)** to ensure complete disruption of the LNPs and full release of encapsulated RNA.

**Note (!):**

- Ensure sufficient mixing during incubation to promote efficient LNP lysis.
- Incomplete lysis may lead to underestimation of total RNA and inaccurate EE% values.

### 4.3. Fluorescence measurement – Total RNA

16. Measure fluorescence again using the same instrument settings as in Step 11.

- ❖ This second measurement corresponds to **total RNA (free + encapsulated)**.

## 5. Data analysis

### 5.1. Background correction

17. Subtract the fluorescence signal of the **blank (no RNA)** from all wells.

### 5.2. Standard curve generation

18. Plot the corrected fluorescence values of the RNA standards versus their concentration. Generate a calibration curve (linear regression).

Notes (!):

- Only use data points within the **linear range** of the assay.
- The **coefficient of determination ( $R^2$ )** should be close to 1 for reliable quantification.




Figure 2. Representative high-range RNA standard curve generated using the RiboGreen assay. Background-corrected fluorescence intensity is plotted as a function of RNA concentration, and a linear regression is applied to determine the calibration curve used for RNA quantification.

### 5.3. Determination of RNA concentrations

19. Use the standard curve to calculate RNA concentrations for:

- Free RNA (before Triton addition)
- Total RNA (after Triton addition)

### 5.4. Encapsulation efficiency (EE%) and yield (EY%) calculation

20. Calculate encapsulation efficiency using Eq.1.

Notes (!):

- EE% is calculated from measurements performed on the **same sample**.
- EE% values may vary depending on the stage at which they are measured (e.g., before or after purification), due to potential changes in LNP integrity or RNA release.

21. (Optional) If RNA input solution was included in the plate, calculate **encapsulation yield** using *Eq.2*.

**Notes (!):**

- **Correct for dilution factors** when comparing RNA input and RNA-LNP samples, especially if different dilutions were applied during plate preparation.
- Take into account the **flow rate ratio (FRR)** used during microfluidic mixing, as it directly affects the final RNA concentration in the RNA-LNP samples.
- EY% may vary depending on the stage at which it is measured (e.g., before or after purification steps such as dialysis or ultrafiltration).

## 5.5. Data processing tools

22. For streamlined analysis, **Inside Therapeutics** provides an **Excel-based calculation tool** that allows:

- EE% and EY% calculations
- Integration of dilution factors and FRR
- Standardized data processing across experiments


## 6. Quality control & troubleshooting

- **Signal range:** Ensure sample fluorescence values fall within the **standard curve range**.
  - If signals are too high → **increase dilution**
  - If signals are too low → **reduce dilution**
- **Reproducibility:** Perform all measurements in **triplicate** and check variability between replicates.
  - High variability may indicate pipetting errors, incomplete mixing, or inconsistent sample preparation.
- **Assay sensitivity:** Maintain consistent **buffer conditions, sample preparation, measurement parameters** across all samples. Small variations in these factors may impact fluorescence intensity and lead to inaccurate quantification.

# Example EE% Data

Representative encapsulation efficiency (EE%) results for standard RNA-LNP formulations prepared using TAMARA microfluidic platform are shown in **Figure 3**. Formulations were produced using CordenPharma LNP Starter Kits (SM-102, ALC-0315, and LP-01) at total flow rates (TFR) of 5 and 10 mL/min, at a fixed flow rate ratio (FRR) of 3. Following formulation, samples were purified by dialysis and analyzed using the RiboGreen assay.

Encapsulation efficiencies were consistently high across all tested conditions, with values generally **around or above 90%** (>80% for all formulations). Comparable EE% values were observed at both flow rates, indicating robust encapsulation performance over the tested processing range.



**Figure 3.** Representative encapsulation efficiency (EE%) of standard RNA-LNP formulations prepared using the TAMARA platform at different total flow rates (TFR = 5 and 10 mL/min). Data are shown for SM-102, ALC-0315, and LP-01 lipid compositions.

## References

1. Quant-iT™ RiboGreen™ RNA Reagent and Kit [Internet]. [cited 2026 Apr 24]. Available from: <https://documents.thermofisher.com/TFS-Assets/LSG/manuals/mp11490.pdf>
2. McKenzie RE, Minnell JJ, Ganley M, Painter GF, Draper SL. mRNA Synthesis and Encapsulation in Ionizable Lipid Nanoparticles. *Curr Protoc*. 2023 Sep 1;3(9). doi:10.1002/cpz1.898 PubMed PMID: 37747354.
3. RNA-LNP Formulation Protocol (SM-102) - Inside Therapeutics [Internet]. [cited 2026 Apr 23]. Available from: <https://insidetx.com/resources/protocols/rna-lnp-formulation-protocol-sm-102/>