

# From Microfluidic Formulation to Payload Quantification: TAMARA and CloudSpec for RNA-LNP Development

RNA-LNP development relies on both **reproducible formulation** and **reliable analytical characterization** to ensure consistent product quality. Encapsulation efficiency (EE%) is a key critical quality attribute, directly influencing dose accuracy, formulation performance, and product consistency.

The **RiboGreen assay** is commonly used to determine EE% but relies on nanoparticle disruption, introducing additional sample preparation and potential variability. **CloudSpec**, based on scatter-free absorbance spectroscopy (SFAS), provides a rapid, lysis-independent alternative for RNA payload quantification.

In this work, SM-102-based polyA-LNPs were formulated using the **TAMARA** microfluidic platform (InsideTx) and characterized using **CloudSpec** (Marama Labs). The effects of the N/P ratio and TFR on particle size, PDI, and EE% were evaluated, while CloudSpec was compared with the conventional RiboGreen assay.



**TAMARA produced highly reproducible RNA-LNPs, while CloudSpec provided rapid and lysis-independent RNA payload quantification.**

*This work was carried out by Aida López Espinar and Kirsty Smith at Marama Labs, in collaboration with Sezen Gul from Inside Therapeutics.*

## Physicochemical characterization of polyA-LNPs

SM-102-based polyA-LNPs were formulated using the TAMARA microfluidic platform under three formulation conditions combining two N/P ratios (6 and 1) and two total flow rates (TFR; 5 and 10 mL/min). Particle size and polydispersity index (PDI) were characterized in triplicate by DLS.

### Key findings

- **LNPs formulated at N/P 6** exhibited highly reproducible particle size (~105 nm) and low PDI (~0.13), independent of TFR (5–10 mL/min).
- **Reducing the N/P ratio from 6 to 1** increased both particle size (~212 nm) and PDI (~0.21), reflecting the reduced availability of ionizable lipid to drive efficient nanoparticle self-assembly.
- **Low variability across replicate batches** confirmed the robustness and reproducibility of the TAMARA microfluidic formulation platform.

TAMARA enables robust and reproducible RNA-LNP formulation while supporting systematic formulation optimization and process development.





Figure 1: Effect of N/P ratio and TFR on particle size and PDI of SM-102-based polyA-LNPs formulated using TAMARA.

# Payload quantification: CloudSpec vs RiboGreen



RNA payload quantification of SM-102 polyA-LNPs was performed using the conventional RiboGreen (RG) assay and the CloudSpec-RiboGreen (CS-RG) assay. Free polyA, total polyA, and EE% were measured for each formulation condition to compare the analytical performance of both methods.

## Key findings

- CloudSpec and the conventional RiboGreen assay showed excellent agreement for free polyA quantification and encapsulation efficiency across all formulation conditions.
- CloudSpec measured total polyA concentrations closer to the expected formulation input, while preserving the same trends between formulations and eliminating the need for nanoparticle lysis.
- Both methods accurately differentiated formulations with high and low encapsulation performance, confirming the effect of decreasing the N/P ratio on RNA encapsulation.




Figure 2: Comparison of RNA payload quantification using the conventional RiboGreen (RG) assay and the CloudSpec-RiboGreen (CS-RG) assay. (a) Free polyA concentration, (b) total polyA concentration, and (c) encapsulation efficiency (EE%) for SM-102-based polyA-LNPs.




Figure 3: CloudSpec measurement principle for lysis-independent RNA payload quantification. The integrating sphere collects transmitted, scattered, and fluorescence-emitted light, enabling determination of total RNA by scatter-free absorbance spectroscopy (SFAS) and free RNA by fluorescence. A representative spectrum and the corresponding EE% calculation are shown.

CloudSpec enables rapid, lysis-independent RNA payload quantification, providing free RNA and EE% measurements consistent with RiboGreen while measuring total RNA closer to the expected formulation input.

## Conclusion

This application note demonstrates the complementary value of TAMARA for reproducible RNA-LNP formulation and CloudSpec for robust RNA payload quantification.

- TAMARA produced highly reproducible RNA-LNPs and enabled systematic evaluation of key formulation parameters, including the N/P ratio and TFR.
- CloudSpec generated EE% measurements consistent with the conventional RiboGreen assay while providing lysis-independent total RNA quantification closer to the expected formulation input.

Together, these complementary technologies support a more efficient and reliable approach to RNA-LNP development, from formulation to payload quantification.



Read the full  
Application Note on  
<https://insidetx.com>